

B6.129S-Tgfbr3^{tm1.1Jba} Genotyping

Oligonucleotides were positioned in exons 2 and 5 to determine the presence or absence of exon 3 across genotypes.

Forward 5'GCTACACCCGACTTGCCCACT3'
Reverse 5'GACCACAGAACCCTCCGAAACC3'

Products were electrophoresed on a 1% agarose/TAE gel and can be seen in Fig. S1D in the Supplementary Data section of Compton, et. al.

Sequencing of the allele is also offered by *Transnetyx*.